

To measure the presence of IgG, IgM, or IgA antibodies against SARS-CoV-2 antigens in serum or plasma samples.

Kit Components

The Guava[®] SARS-CoV-2 Multi-Antigen Antibody Assay (FCPA100101) includes the following components:

- Guava SARS-CoV-2 Microsphere Mix (14-20327, 5.5 mL)
- PBS-TBN Wash Buffer (14-20335, 125 mL)
- Anti-Human IgG Detection Reagent (14-20332, 5.5 mL)
- Anti-Human IgM Detection Reagent (14-20334, 5.5 mL)
- Anti-Human IgA Detection Reagent (14-20333, 5.5 mL)

Storage Conditions

- Store Guava[®] SARS-CoV-2 Multi-Antigen Antibody Kit refrigerated at 2°C to 8°C.
- Do not freeze.
- Protect the Guava SARS-CoV-2 Multi-Antigen Antibody Assay from excessive exposure to light.

Recommended Materials

- Guava[®] Muse[®] Cell Analyzer, Guava[®] easyCyte™ flow cytometer, or other flow cytometer equipped with a 488 nm or 532 nm laser.
- Magnetic base rack or base plate for manual washing.
- Micropipettors
- Disposable micropipettor tips
- Vortex mixer with adaptor
- Orbital plate shaker
- Timer
- Deionized (DI) water
- 10% bleach solution
NOTE: 10% bleach is defined as 0.6% sodium hypochlorite
- 1.5 mL sample tubes or 96-well plate

Assay Protocol

- Dilute serum or plasma samples 1:400
- ↓
- Add 50 µL of control or diluted sample to each tube or well.
- ↓
- Add 50 µL of Guava[®] SARS-CoV-2 Microsphere Mix to each tube or well.
- ↓
- Seal the plate or cap the tube and incubate while shaking at room temperature for 30 minutes.
- ↓
- Remove the liquid and wash the 1X well on the magnet with 150 µL of PBS-TBN Wash Buffer.
- ↓
- Add 50 µL of Anti-Human Detection Reagent.
- ↓
- Seal the plate or cap the tube and incubate while shaking at room temperature for 30 minutes.
- ↓
- Remove the liquid on the magnet.
- ↓
- Resuspend the microspheres in 200 µL of PBS-TBN Wash Buffer for reading on the flow cytometer.



Guava[®] SARS-CoV-2 Multi-Antigen Antibody Assay Quick Reference Card

Expected Results

The following figures show examples of typical results obtained using the Guava[®] SARS-CoV-2 Antibody Assay on the Guava[®] Muse[®] Cell Analyzer (Figure 1) and with the Guava[®] easyCyte[™] flow cytometer (Figure 2).




Figure 1: Results from a SARS-CoV-2 positive human serum sample. The sample was prepared with the Guava SARS-CoV-2 Multi-Antigen Antibody Assay and acquired on the Guava Muse Cell Analyzer using the SARS-CoV-2 Antibody module. The statistics show the count of the bead population and median fluorescent intensity for each population.




Figure 2: Results from a SARS-CoV-2 positive human serum sample. The sample was prepared with the Guava SARS-CoV-2 Multi-Antigen Antibody Assay and acquired on the easyCyte 12HT flow cytometer using InCyte[™] software. Output for statistics of bead counts and median fluorescent intensity of bead populations can be set up and obtained using InCyte software.

This document is a supplemental tool and consistent with the instructions for use in the Guava[®] SARS-CoV-2 Multi-Antigen Assay Package Insert (89-30000-00-959). Please refer to the package insert for additional information. For more information on Guava flow reagents and kits, refer to <https://www.luminexcorp.com/flow-cytometry-kits-and-reagents/>.