



complexity simplified.

Quick Reference Card

Muse[®] Human CD8 T Cell Kit MIM100102

For use in laboratory research studies to determine the CD8 T cell count, total lymphocyte count, and CD8 cell percent of lymphocytes in human whole blood and PBMC samples

CAUTION: This kit should not be used for patient diagnosis or patient management.

For Research Use Only. Not for use in diagnostic procedures.

Storage Conditions

Store the Muse[®] Human CD8 T Cell Cocktail and 1X Assay Buffer BA at 2 to 8°C, protected from light. Store the Human 1X Lysing Solution at room temperature.


Kit Components

- Muse[®] Human CD8 T Cell Cocktail (Part No. 4700-1625, 100 tests/vial)
- Human 1X Lysing Solution (Part No. 4700-1620, 100 tests/vial)
- 1X Assay Buffer BA (Part No. 4700-1360, 100 tests/vial)

Materials Recommended

- Guava[®] Muse[®] Cell Analyzer
- Whole blood or PBMC samples
- Micropipettors
- Disposable micropipettor tips
- Microcentrifuge tubes with screw caps, 1.5 mL (VWR Catalog No. 16466-030, or equivalent)
- Vortex mixer

Assay Protocol



* PBMC samples should be in 1X Assay Buffer BA at a concentration of 4×10^6 to 2×10^7 cells/mL.

NOTE: A detailed kit user's guide can be found at www.luminexcorp.com/flowkits.


Expected Results

The Muse® Human CD8 T Cell Assay provides results for CD8 T lymphocyte count, CD8 percent of total lymphocytes, and total lymphocyte count.

Results for each of the four quadrants are as follows:

- Upper-left quadrant: CD8(-) T&B lymphocytes [T&B lymphocyte(+) and CD8(-)]
- Upper-right quadrant: CD8(+) T lymphocytes [T&B lymphocyte(+) and CD8(+)]
- Lower-left quadrant: negative cells [T&B lymphocyte(-) and CD8(-)]
- Lower-right quadrant: CD8(+) cells [T&B lymphocyte(-) and CD8(+)]

Figures A and B: Example Data: Results obtained with the Muse Human CD8 T Cell module for adult whole blood stained with the Muse Human CD8 T Cell Kit and acquired on the Guava® Muse Cell Analyzer. Figure A shows results without dot plots, while Figure B shows results with optional dot plots. The statistics show cells/ μL in the original sample, population percentages, as well as total lymphocyte count and CD8 T cell percentage of lymphocytes. The first plot in Figure B shows Cell Size vs Lymphocytes and the second plot shows CD8 staining vs Lymphocytes.



The latest version of Muse® software, which includes all assay modules, as well as the kit user's guide, can be found at www.luminexcorp.com/flowkits.

Related Products

For Research Use Only. Not for use in diagnostic procedures.

Muse® Human CD4 T Cell Kit – MIM100101

Muse® Human B Cell Kit – MIM100103

Muse® Human Lymphocyte CD69 Kit – MIM100104

Muse® Human Lymphocyte CD25 Kit – MIM100105

Muse® System Check Kit – MCH100101

Muse® Count & Viability Kit (40mL) – MCH100102

Muse® Annexin V & Dead Cell Kit – MCH100105

Muse and Guava are trademarks of Luminex Corporation, registered in the U.S. and other countries..

Part No. 4600-3438 Rev D, February 2020. Printed in the USA.

© 2020 Luminex Corporation. All rights reserved.