



Cytek® cFluor® Human Pan Leukocyte Kit, LNW

Cytek® cFluor® Human Pan Leukocyte Kit, LNW helps researchers to identify and enumerate leukocyte subsets in drug discovery and development, and their other research needs. This has been designed to fully enumerate the complete set of major leukocyte subsets and to mirror and expand on those identified in a traditional complete white blood cell differential. The kit has been optimized to be run on whole blood as a lyse no wash (LNW) assay which, in tandem with a complete blood count (CBC), allows accurate dual platform absolute population counts. Furthermore, the kit has been designed and validated to allow the sequential identification of populations that provides for easy and highly reproducible gating of results. This kit can also serve as a backbone panel which can be expanded upon to further identify any populations of interest, as well as the functional or activation/exhaustion status of a particular population.

PRODUCT DETAILS																																																	
Catalog number:	R7-40007 (50 Tests)																																																
Category:	Immunophenotyping																																																
Format:	cFluor® conjugated antibodies in individual vials																																																
	<table border="1"> <thead> <tr> <th>Target</th> <th>Clone</th> <th>Fluorochrome</th> </tr> </thead> <tbody> <tr><td>CD8</td><td>SK1</td><td>cFluor® V450</td></tr> <tr><td>HLA-DR</td><td>L243</td><td>cFluor® V505</td></tr> <tr><td>CD45</td><td>HI30</td><td>cFluor® V547</td></tr> <tr><td>CD4</td><td>SK3</td><td>cFluor® V610</td></tr> <tr><td>CD16</td><td>3G8</td><td>cFluor® B515</td></tr> <tr><td>CD34</td><td>4H11</td><td>cFluor® BYG575</td></tr> <tr><td>CD123</td><td>6H6</td><td>cFluor® BYG610</td></tr> <tr><td>CD193 (CCR3)</td><td>5E8</td><td>cFluor® BYG667</td></tr> <tr><td>CD56</td><td>LT56</td><td>cFluor® BYG710</td></tr> <tr><td>CD19</td><td>HIB19</td><td>cFluor® BYG750</td></tr> <tr><td>CD14</td><td>MEM-15</td><td>cFluor® BYG781</td></tr> <tr><td>CD7</td><td>CD7-6B7</td><td>cFluor® R659</td></tr> <tr><td>CD20</td><td>2H7</td><td>cFluor® R685</td></tr> <tr><td>CD66b</td><td>G10F5</td><td>cFluor® R720</td></tr> <tr><td>CD3</td><td>SK7</td><td>cFluor® R780</td></tr> </tbody> </table>	Target	Clone	Fluorochrome	CD8	SK1	cFluor® V450	HLA-DR	L243	cFluor® V505	CD45	HI30	cFluor® V547	CD4	SK3	cFluor® V610	CD16	3G8	cFluor® B515	CD34	4H11	cFluor® BYG575	CD123	6H6	cFluor® BYG610	CD193 (CCR3)	5E8	cFluor® BYG667	CD56	LT56	cFluor® BYG710	CD19	HIB19	cFluor® BYG750	CD14	MEM-15	cFluor® BYG781	CD7	CD7-6B7	cFluor® R659	CD20	2H7	cFluor® R685	CD66b	G10F5	cFluor® R720	CD3	SK7	cFluor® R780
Target	Clone	Fluorochrome																																															
CD8	SK1	cFluor® V450																																															
HLA-DR	L243	cFluor® V505																																															
CD45	HI30	cFluor® V547																																															
CD4	SK3	cFluor® V610																																															
CD16	3G8	cFluor® B515																																															
CD34	4H11	cFluor® BYG575																																															
CD123	6H6	cFluor® BYG610																																															
CD193 (CCR3)	5E8	cFluor® BYG667																																															
CD56	LT56	cFluor® BYG710																																															
CD19	HIB19	cFluor® BYG750																																															
CD14	MEM-15	cFluor® BYG781																																															
CD7	CD7-6B7	cFluor® R659																																															
CD20	2H7	cFluor® R685																																															
CD66b	G10F5	cFluor® R720																																															
CD3	SK7	cFluor® R780																																															
Test Dilution:	2 µl per test																																																
Application:	Flow cytometry																																																
Formulation:	Phosphate-buffered saline, pH 7.2, containing 0.09% sodium azide and 0.2% BSA (BSA Country of Origin USA)																																																
Storage:	2-8°C and protected from light. Do not freeze																																																

PRODUCT DESCRIPTION

The Cytex® cFluor® Human Pan Leukocyte Kit, LNW allows for the identification of more than 20 leukocyte subpopulations, including neutrophils, eosinophils, basophils, lymphocytes, hematopoietic stem cells, monocyte subsets, T cell subsets, B cells, and NK cell subsets in whole blood.

CD3 is expressed on all mature T cells, NK T cells, and some thymocytes. CD3, a part of the CD3/T cell receptor complex, plays a role in antigen recognition, signal transduction, and T cell activation.

CD4 is expressed on most thymocytes, a major subset of T cells, and on monocytes/macrophages. Functionally, CD4 is associated with thymic differentiation, in conjunction with MHC class II molecules in antigen recognition and with signal transduction.

CD7, also known as GP40, is found on T cells, NK cells, thymocytes, and pre-B cells. CD7 is called T cell leukemia antigen for its expression on acute lymphocytic leukemia and some acute myeloid leukemia.

R8-50007 Rev. A, Effective Date: 01/13/2023

Cytex Biosciences, Inc. | 47215 Lakeview Boulevard, Fremont, CA 94538

Phone: +1 (510) 657-0102 | Fax: +1 (510) 657-0151 | cytekbio@cytekbio.com | www.cytekbio.com

CD8 is found on thymocytes, on a subset of T cells and on NK cells. This molecule acts as a co-receptor in MHC class I molecules in antigen recognition and plays a role in T cell activation and thymic differentiation.

CD14 is highly expressed on monocytes and macrophages and is known as a high affinity LPS receptor. It is also expressed on granulocytes, but at a lower level. In addition, CD14 is found on interfollicular dendritic cells, reticular dendritic cells, and Langerhans cells.

CD16 is expressed on NK cells, monocytes and macrophages in the form of CD16a. Another form of CD16, CD16b is expressed on neutrophils. CD16 engagement of IgG leads to NK cell activation, antibody-dependent cell-mediated cytotoxicity (ADCC) and phagocytosis

CD19 is expressed in the B cell lineage, from pro-B to blastoid B cells. However, it is absent on plasma cells. It is also expressed on follicular dendritic cells. CD19 is involved in B cell development, activation, and differentiation.

CD20 is a membrane protein specific to mature B cells and is involved in the differentiation of B cells into plasma cells. The expression of CD20 is also found on a variety of malignant B cells.

CD34, also known as mucosialin, is a transmembrane protein that is commonly expressed on hematopoietic stem cells. The CD34 is also expressed on some populations of mesenchymal stem cells and vascular endothelium. The stem cell antigen CD34 is functioning as a regulator of hemopoietic cell adhesion.

CD45 is expressed on all hematopoietic cells, except erythrocytes and platelets. CD45 is a signaling molecule that is involved in cellular proliferation, differentiation and in regulation of immune cell functions.

CD56, within the hematopoietic system is expressed on NK cells and NKT cells, a subset of T cells. In the nervous system, CD56 is expressed by neurons and plays a role in the homotypic adhesion of neural cells.

CD66b, also known as CD67, CGM6, and NCA-95, is exclusively expressed on human neutrophils and eosinophils. It is also recognized as granulocyte "activation marker".

CD123, also known as the alpha chain of the IL-3 receptor, is expressed predominantly on basophils, dendritic cells, hematopoietic progenitor cells, macrophages, eosinophils, monocytes, megakaryocytes, and some B cells.

HLA-DR is present on the surface of antigen-presenting cells, including B cells, dendritic cells, macrophages, monocytes and activated T cells. MHC class II regulates the immune system by playing a critical role in binding and presenting antigen-derived peptides to peptide-MHC class II-specific CD4 T cells.

CD193 (CCR3) is expressed on eosinophils, basophils, mast cells, and T cells. CCR3 is a major receptor for eotaxin, RANTES, and MCP-4. These chemo-attractants are involved in allergic inflammation and infiltration of CCR3 expressing cells into inflamed tissues.

REFERENCES

- Clevers H, et al. Ann Rev Immunol. 1988 **6**, 629 (1988)
Evans RL, et al. Immunol. **78**, 544 (1981)
Rabinowich H. J Immunol. **152**(2):517-526 (1994)
Cole DK, et al. Cell Mol Immunol. **1**, 81 (2004)
Ferrero E, et al. J Immunol. **145**, 331 (1990)
Wirthmueller U, et al. J Exp Med. **175**, (1992)
Passlick B, et al. 1989. Blood. **74**(7):2527-34.
Macardle PJ, et al. Cell Immunol. **168**, 59 (1996)
Pieper K, et al. J Allergy Clin Immunol. **131**, 959-971 (2013)
Berenson RJ, et al. Blood. **77**(8):1717-22. (1991)
Healy L, et al. PNAS. 92(26):12240-12244. (1995)
Loken MR, et al. Cytometry. **11**, 453 (1990)
Poli A, et al. Immunology. **126**, 458 (2009)
Yoon J, et al. J Immunol. **179** (12) 8454-8462 (2007)
Yamada T, et al. Allergy Clin Immunol. **101**, 677 (1998)
Uguccioni M, et al. J Clin Invest. **1**;100(5):1137-43 (1997)
Walseng E, et al. J Biol Chem. **283**,14717 (2008)

RECOMMENDED USAGE

Whole blood collected in K₂EDTA, Heparin, ACD and Cyto-Chex[®] BCT blood tubes have been tested to validate the performance of this kit. For staining procedures, product data, and gating strategy, please refer to the Reagents and Protocols sections of our website at www.cytokbio.com.

Please briefly centrifuge the reagent vial before use.

Use appropriate personal protective equipment per the product safety data sheet when using this product.



For Research Use Only. Not intended for use in diagnostic procedures.

cFluor® V547, cFluor® B515, cFluor® R685 and cFluor® R720 are equivalent to CF®405L, CF®488A, CF®660C and CF®700 respectively, manufactured and provided by Biotium, Inc. under an Agreement between Biotium and Cytek (LICENSEE). The manufacture, use, sale, offer for sale, or import of the product is covered by one or more of the patents or pending applications owned or licensed by Biotium. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel.

cFluor® BYG610, cFluor® BYG667, cFluor® BYG710, cFluor® BYG750, and cFluor® BYG781 are a tandem dyes made with R-PE. cFluor® R780 is a tandem dye made with APC. Caution – Tandem dyes may show changes in their emission spectra with prolonged exposure to light or fixatives.

"Cytek" and "cFluor" are trademarks or registered trademarks of Cytek Biosciences, Inc. All other service marks, trademarks and tradenames appearing herein are the property of their respective owners.