

Acquisition Protocol for cFluor™ 14 Color Immunoprofiling Kit


This acquisition protocol provides step-by-step instructions to set up your Cytek® Aurora or Northern Lights system (3-laser V-B-R configuration or higher) for data acquisition of the cFluor™ 14 Color Immunoprofiling Kit (R7-40000). This protocol provides instructions on: 1) preparing the SpectroFlo® software; 2) setting up the instrument; 3) acquiring controls and samples. This kit contains 25 tests.

NOTE: For suggestions on how to prepare human PBMC or whole blood samples, download the protocol from the cFluor™ kit [product page](#).

Preparing the SpectroFlo® Software Adding the 14 cFluor dyes to the library

NOTE: Adding cFluor dyes to the library is only required for SpectroFlo® versions 2.2 and below.

1. Download the cFluor fluorescent tag list CSV from the cFluor™ kit [product page](#).
2. Import this CSV file into the **Library** module, under **Fluorescent Tags** to add in the 14 cFluor™ fluorochromes used in this kit.




The screenshot shows the SpectroFlo software interface. The 'Library' module is active, and the 'Fluorescent Tag Groups' table is displayed. The 'Import' button in the top toolbar is highlighted with a red box. The 'Fluorescent Tags' menu item in the left sidebar is also highlighted with a red box. The table contains the following data:

Name	Created By	Date Created	Name	Emission Wavelength	Laser Excitation	Display Name	Synonyms
Blue Laser	System	December 04, 2018 - 16:06 PM	cFluor V420	420	405	cFluor V420	cFluorV420
Red Laser	System	December 04, 2018 - 16:06 PM	cFluor V450	450	405	cFluor V450	cFluorV450
Violet Laser	System	December 04, 2018 - 16:06 PM	cFluor B515	515	488	cFluor B515	cFluorB515
UV Laser	System	December 04, 2018 - 16:06 PM	cFluor V547	547	405	cFluor V547	cFluorV547
YellowGreen Laser	System	December 04, 2018 - 16:06 PM	cFluor BYG575	575	561	cFluor BYG575	cFluorBYG575
Fluorescent Proteins	System	December 04, 2018 - 16:06 PM	cFluor R659	659	635	cFluor R659	cFluorR659
Viability	System	December 04, 2018 - 16:06 PM	cFluor BYG667	667	561	cFluor BYG667	cFluorBYG667
New dyes	Admin	March 11, 2019 - 14:14 PM	cFluor R668	668	635	cFluor R668	cFluorR668
cFluor Immunoprofiling Kit Tags	Admin	February 22, 2021 - 16:23 PM	cFluor B690	690	488	cFluor B690	cFluorB690
			cFluor BYG710	710	561	cFluor BYG710	cFluorBYG710
			cFluor R720	720	635	cFluor R720	cFluorR720
			cFluor R780	780	635	cFluor R780	cFluorR780
			cFluor BYG781	781	561	cFluor BYG781	cFluorBYG781
			cFluor R840	840	635	cFluor R840	cFluorR840

NOTE: If any cFluors were previously entered into the SpectroFlo® Library, a warning message will appear. Click **Overwrite** to overwrite the information in the library with the new fluorescent tag information.

Import the experiment template

1. Download the cFluor Immunoprofiling Kit experiment template from the kit [product page](#). This template includes a reference group with predefined stopping criteria, assigned marker names, as well as recommended acquisition and analysis worksheets.
2. Import the template into the **Library** module, under **Experiment Templates**.




The screenshot shows the 'Library' interface with the 'Experiment Templates' section selected. The 'Library' tab in the top navigation bar is highlighted with an orange box. The 'Experiment Templates' section in the left sidebar is also highlighted with an orange box. The main content area shows a table with one entry: 'Default'.

Name	Created By	Date Created	Description
✓ Default	Admin	April 08, 2020 - 11:49 AM	Default Experiment

Setting up the Instrument

1. Follow the instructions for instrument startup and **Performing Daily QC** as outlined in the User's Guide.
2. From the **Acquisition** module, create a new experiment by clicking on **Template**, then choosing "14C Immunophenotyping Kit".



3. For best results, we recommend running cells for single color controls. If desired, control type can be changed to beads.
4. The sample tube is set to acquire 500,000 lymphocytes. This can be changed by clicking **Edit**, then changing the Stopping Gate and Stopping Criteria under **Acquisition**.
5. Add sample tubes and groups as needed. Marker labels will auto-populate for any tubes added to the existing sample group. To preserve these markers for additional groups, duplicate the existing group and rename the tube as necessary.
6. Preview unstained cell control at low flow rate to minimize wasted sample volume. Starting from the default CytekAssaySettings, optimize the FSC and SSC gains, as well as the FSC threshold (see Figure 1).

NOTE: Instrument settings can be saved as "cFluor IP Kit - 14C" for future use by clicking the **Save As** button in the **Instrument Control** window.






Figure 1: Example cell resolution of human PBMC samples after FSC and SSC adjustment.

Acquiring Controls and Samples


1. When running unstained and single color controls, acquire all tubes in the **Reference Control** group. Click **Unmix**. **NOTE:** Refer to Appendix A for additional workflows when storing and reusing single color controls.



2. For each single color control, gate on the appropriate population in the scatter plot and adjust the positive and negative gates in the histograms. Ensure that the signatures for the single color controls match the expected spectrum. Refer to Appendix B for the expected spectra of each cFluor and placement of positive and negative gates in the histograms.
3. Click **Next**. Under the **QC Controls** tab, click on **Similarity Matrix** to visualize how similar the 14 cFluors are in the panel. Click on **View Similarity Index** to compare to expected values found in Appendix C.



4. Click **Live Unmixing**.
5. Acquire all multicolor samples. Adjust the gates in the analysis worksheet as needed.
6. To export the FCS files for offline analysis: 1) close the experiment; 2) click **My Experiments**; 3) select the experiment and click **Export**; 4) choose a directory and click **Export**.



Appendix A: Reusing Single Color Controls


To store controls for reuse in future experiments, follow the instructions below.


NOTE: For best results, the instrument should be maintained properly, QC performed daily and the same lot of reagents used.


Reusing Single Color Controls From a Previous Experiment

1. From the **Acquisition** module, open **My Experiments**. Right click the saved experiment and select **Duplicate**. This will duplicate the data.
2. Keep all of the FCS files in the **Reference Control** group.
3. Add, delete, duplicate, and/or rename multicolor groups and samples as necessary.
NOTE: Marker labels may need to be added for new sample groups and tubes.
4. Acquire all multicolor samples. Adjust the gates in the analysis worksheet as needed. Unmixing will carry over from the previous experiment when references were initially run.
5. To export the FCS files for analysis at a different workstation: 1) close the experiment; 2) click **My Experiments**; 3) select the experiment and click **Export**; 4) choose a directory and click **Export**.

Appendix B: Single Color Control Signatures for 3-Laser (V-B-R) Cytek Aurora and Northern Lights and Suggested Gating When Using Cells as Controls







NOTE: For single color control signatures for other Cytek Aurora and Northern Lights configurations please refer to the fluorochrome guidelines on the [Resources Page](#).

NOTE: For cFluor V450 and cFluor V547, the Experiment Template has been formatted to assign the Unstained Cells as the Negative Control, since no negative populations are present when gating on monocytes (CD14 cFluor V450) and lymphocytes (CD45 cFluor V547) for these two controls.

