




# cFluor<sup>®</sup> Anti-Human CD45

## Instructions For Use

Catalog No.	Test/Vial	Product Name
R7-11009	100	cFluor <sup>®</sup> B690 Anti-Human CD45 (2D1)
R7-11010	25	cFluor <sup>®</sup> B690 Anti-Human CD45 (2D1)
R7-11011	100	cFluor <sup>®</sup> V547 Anti-Human CD45 (HI30)
R7-11012	25	cFluor <sup>®</sup> V547 Anti-Human CD45 (HI30)

## Copyright and Trademarks

© 2022 Cytex Biosciences, Inc. All rights reserved. Cytex, the Cytex logo, cFluor and Northern Lights are trademarks or registered trademarks of Cytex Biosciences, Inc. All other service marks, trademarks and tradenames are the property of their respective owners.

 Cytex Biosciences, Inc  
47215 Lakeview Blvd.  
Fremont, CA 94538  
USA  
1.877.92.CYTEK (1.877.922.9835)

products@cytekbio.com  
cytekbio.com



Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague  
Netherlands

## 1. Intended use

This product is intended for in vitro diagnostic use to identify human cells expressing CD45 antigen molecules in countries where the regulatory approval has been obtained from the local regulatory authorities.

## 2. Application

The CD45 monoclonal antibody binds to isoforms of human CD45, a 180-240-kDa type I membrane glycoprotein, also known as LCA (leukocyte common antigen). CD45 is the first and prototypic receptor-like protein tyrosine phosphatase that expresses on all human leukocytes. It is absent on mature erythrocytes, platelets, and non-hematopoietic cells. The antibody is conjugated to a fluorophore and purified by affinity chromatography.

## 3. Components

CD45 monoclonal antibody conjugated with the following listed cFluor fluorescent dye is supplied in phosphate-buffered saline, pH 7.2, containing 0.09% sodium azide and 0.2% BSA (BSA Country of Origin USA).

Antibody specificity	CD45	CD45
Clone	2D1	HI30
Immunoglobulin subtype	IgG1, kappa	IgG1, kappa
Species and genus	Mouse	Mouse
Fluorescent dye	cFluor <sup>®</sup> B690 <sup>2</sup>	cFluor <sup>®</sup> V547 <sup>1</sup>
Excitation wavelength	488 nm	405 nm
Emission peak	690 nm	547 nm

## 4. Storage and Handling

This product is stable until the expiration date shown on the label when stored away from light at 2 ~ 8 °C. Do not freeze.

## 5. Other Materials required but not supplied

- RBC lysing solution
- Pipettes and pipette tips of 20 µL, 100 µL and 1000 µL
- 12x75mm tube
- Vortex mixer
- Flow cytometer

## 6. Specimen Requirements

- 1 Require peripheral blood of not less than 500 µL collected by venipuncture in EDTA anticoagulation tube.
- 2 After collection, the samples should be stored at room temperature (18 ~ 25 °C). Avoid shaking. The storage time should not exceed 24 hours.

- 3 After staining, the samples should be stored at 2 ~ 8 °C away from light and analyzed by flow cytometry within 2 hours.
- 4 Avoid samples with microbial contamination or coagulation.

## 7. Procedure

- 1 Add 100 µL well-mixed EDTA anticoagulated whole blood to the bottom of a tube. Avoid blood touching the upper tube wall.
- 2 Briefly centrifuge this product before use. Add 5 µL of CD45-cFluor-conjugated reagent to the bottom of the tube.
- 3 Mix well by vortex and incubate for 15-30 minutes at room temperature and away from light.
- 4 Add 2 mL of 1 X lysis buffer into the tube, mix briefly by vortex, and incubate for 10-15 minutes at room temperature in the dark.
- 5 Centrifuge at 300g for 5 minutes, discard the supernatant, add 2 mL PBS with 0.02% BSA, and 0.09% NaN<sub>3</sub> to resuspend the cell.
- 6 Centrifuge at 300g for 5 minutes, discard the supernatant, add 300 µL PBS with 0.02% BSA, and 0.09% NaN<sub>3</sub> to resuspend the cells and keep at 4 °C, and analyze on flow cytometer within 2 hours. If delayed analysis is needed (more than 2 hours), 300 µL of PBS containing 1% paraformaldehyde should be used to resuspend the cells and store the sample in a refrigerator at 2-8 °C away from light, but the storage time should not exceed 24 hours.

## 8. Quality Control

- Instrument QC: Use the manufacturer recommended controls according to the model of the flow cytometer.
- Refer to the instrument User's Guide for instrument maintenance.

## 9. Warnings

- This reagent contains traces of sodium azide. Do not pipette by mouth.
- Use appropriate personal protective equipment per the safety data sheet when using this product.
- Follow biosafety practice in compliance with federal, state, and local regulations to handle all biological samples and materials in contact with them.
- Contact Cytex Support or refer to [cytekbio.com](http://cytekbio.com) for details on troubleshooting.

## 10. Performance Characteristics

### 10.1. Accuracy

Three replicate tubes were stained with each CD45-cFluor-conjugated reagent and analyzed on Cytek Northern Lights™ flow cytometer. The percent CD45+ lymphocytes results were within the control blood target value range provided by the manufacturer.

Specimen: CD-CHEX PLUS	Percent CD45+ Lymphocytes				
CD45-Fluorescent Dye	R1	R2	R3	Mean	Target Value Range
cFluor B690	32.4	33.4	31.6	32.5	22.0-42.0
cFluor V547	31.1	31.2	31.2	31.2	20.9-40.9

### 10.2. Intra-batch precision

Ten replicate tubes were stained with the same batch of each CD45-cFluor-conjugated reagent and analyzed on Cytek Northern Lights™ flow cytometer. The CV of percent CD45+ lymphocytes was calculated and was within the acceptance criteria.

Specimen: Normal Blood	Percent CD45+ Lymphocytes		
CD45-Fluorescent Dye	Average (%)	% CV	Criteria
cFluor B690	31.7	1.18	CV≤8%
cFluor V547	18.5	4.19	CV≤15%

### 10.3. Inter-batch precision

Three replicate tubes were stained with three batches of each CD45-cFluor-conjugated reagent and analyzed on Cytek Northern Lights™ flow cytometer. The CV of percent CD45+ lymphocytes was calculated and was within the acceptance criteria.

Specimen: CD-CHEX PLUS	Percent CD45+ Lymphocytes		
CD45-Fluorescent Dye	Average (%)	% CV	Criteria
cFluor B690	31.8	1.97	CV≤8%
cFluor V547	30.7	4.52	

## 10.4. Staining stability

Three replicate tubes were stained with the same batch of each CD45-cFluor-conjugated reagent and analyzed on Cytek Northern Lights™ flow cytometer at these timepoints: within 2-hour (T0), 6-hour, 24-hour, 48-hour, and 72-hour (for V457 only) after staining. The percent CD45+ lymphocytes at each time point were compared to T0, and the mean relative difference was calculated and was within the acceptance criteria.

Specimen: Normal Blood	Percent CD45+ Lymphocytes					
CD45-Fluorescent Dye	Average (%)	Relative Difference vs. 2H				Criteria
		6H	24H	48H	72H	
cFluor B690	31.5	0.90%	0.47%	-0.06%	NA	Relative Difference ≤10%
cFluor V547	30.3	-0.35%	-0.60%	-1.12%	-5.02%	

## 10.5. Dilution linearity

The samples were serial diluted into five levels (undiluted, 2X, 4X, 8X, 16X). Four replicate tubes at each dilution level were stained with the same batch of each CD45-cFluor-conjugated reagent and analyzed on Cytek Northern Lights™ flow cytometer. The median of percent CD45+ lymphocytes at each dilution level were compared to the median of percent CD45+ lymphocytes at all levels, the relative difference was calculated and was within the acceptance criteria.

Specimen: CD-CHEX-PLUS	Percent CD45+ Lymphocytes					
CD45-Fluorescent Dye	Relative Difference					Criteria
	Undiluted	2X dilution	4 X dilution	8 X dilution	16 X dilution	
cFluor B690	0.08%	1.10%	0.06%	0.06%	-0.67%	Relative Difference ≤10%
cFluor V547	2.63%	2.39%	1.96%	-1.64%	-1.38%	

## 11. Limitations

- 1 This reagent can be used with a flow cytometer and is not recommended for fluorescence microscopy and immunohistochemistry.
- 2 This reagent is a fluorescent labeled product. It is easy to quench with extended light exposure and should be handled away from light.
- 3 If not following the lyse wash procedure described above, the reagent performance can be affected.
- 4 The results may be affected by improper storage of reagents, coagulation of specimens, improper storage of specimens, or incomplete lysis of red blood cells in the samples.
- 5 The test results of this reagent are for clinical reference only. Patient history, other laboratory tests, and treatment response should also be considered for diagnosis.

## 12. References

- Terry LA, et al. 1990. Immunology. 64:331
- Hermiston ML, et al. 2003. Annu Rev Immunol. 21:107

<sup>1</sup>cFluor® V547 is equivalent to CF®405L manufactured and provided by Biotium, Inc. under an Agreement between Biotium and Cytek (LICENSEE). The manufacture, use, sale, offer for sale, or import of the product is covered by one or more of the patents or pending applications owned or licensed by Biotium. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel.

<sup>2</sup>cFluor® B690 is a tandem dye made with PerCP. Caution: Tandem dyes may show changes in their emission spectra with prolonged exposure to light or fixatives.