

cFluor® R720 Anti-Human CD38 (HIT2)

PRODUCT DETAILS	
Catalog Number:	R7-20061 (100 tests) R7-20062 (25 tests)
Lot Number:	F-xxxxxx-xx
Reactivity:	Human, Chimpanzee, Horse
Clone:	HIT2
Format:	cFluor® R720*
Isotype:	Mouse IgG1, κ
Concentration:	50 µg / mL
Application:	Flow cytometry
Formulation:	Phosphate-buffered saline, pH 7.2, containing 0.09% sodium azide and 0.2% BSA (BSA Country of Origin USA)
Storage:	2-8°C and protected from light. Do not freeze
Expiration Date:	Month DD, YYYY

This lot has been tested according to Cytek QC testing protocols and has met Cytek QC test specifications.

Reviewed by

Date

For Research Use Only. Not intended for use in diagnostic procedures.

*cFluor® R720 is equivalent to CF® 700 manufactured and provided by Biotium, Inc. under an Agreement between Biotium and Cytek (LICENSEE). The manufacture, use, sale, offer for sale, or import of the product is covered by one or more of the patents or pending applications owned or licensed by Biotium. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel.

R8-80061 Rev. B, Effective Date: 05/27/2022

Cytek Biosciences Inc. | 47215 Lakeview Boulevard, Fremont, CA 94538

Phone: +1 (510) 657-0102 | Fax: +1 (510) 657-0151 | cytekbio@cytekbio.com | www.cytekbio.com