



## cFluor<sup>®</sup> BYG750 Anti-Human CD4 (SK3)

PRODUCT DETAILS	
<b>Catalog Number:</b>	R7-20159 (100 tests) R7-20160 (25 tests)
<b>Lot Number:</b>	F-xxxxxx-xx
<b>Reactivity:</b>	Human
<b>Clone:</b>	SK3
<b>Format:</b>	cFluor <sup>®</sup> BYG750
<b>Isotype:</b>	Mouse IgG1, $\kappa$
<b>Concentration:</b>	20 $\mu$ g / mL
<b>Application:</b>	Flow cytometry
<b>Formulation:</b>	Phosphate-buffered saline, pH 7.2, containing 0.09% sodium azide and 0.2% BSA (BSA Country of Origin USA)
<b>Storage:</b>	2-8°C and protected from light. <b>Do not freeze</b>
<b>Expiration Date:</b>	Month DD, YYYY

This lot has been tested according to Cytex QC testing protocols and has met Cytex QC test specifications.

\_\_\_\_\_  
Reviewed by

\_\_\_\_\_  
Date

**For Research Use Only. Not intended for use in diagnostic procedures.**